



STARMed Ventumask

For CPAP therapy with Venturi flow driver and adjustable PEEP valve



Now available

with nebuliser port option

CV205S (Small)

CV205M/L (Medium/Large)

CV205XL (Extra Large)

Nebuliser kit not included

Critical Care ▪ Patient Interfaces ▪ **STARMed**


StarMed Ventumask with Venturi flow driver & adjustable PEEP valve

StarMed Ventumask

StarMed Ventumask is a unique design for CPAP therapy that is simply activated by connecting to an appropriate oxygen supply. It can be used in a hospital environment or pre-hospital use, for example in an ambulance.


Advantages and Technical Features

- Ventumask is a complete and easy to use system, where the oxygen supply tubes, a Venturi flow driver and a PEEP valve are integral to the CPAP mask.
- The adjustable PEEP valve is designed exclusively for Ventumask and the administration to the patient is facilitated at some pre-settings: 5 - 7.5 - 10 - 12.5 cmH₂O
- The integrated Venturi flow driver can generate high patient flows while consuming a reduced amount of compressed oxygen
- O₂% can be adjusted from 40% to 100% at each PEEP level
- The user is guided by a performance table integral to the product, for independent PEEP and O₂% settings



The kit is composed of:

- Venturi flow driver integral to the mask
- Adjustable headgear included
- CPAP mask with safety valve and port to monitor actual O₂% or pressure
- Adjustable PEEP valve designed for Ventumask facilitates some pre-calibrated pressure settings
- Integral oxygen tubes



Dual flow meter
Model: STS.15/30
Dual oxygen flow meter
15 + 30 L/min with O₂ coupling

Ventumask for CPAP therapy with Venturi flow driver and adjustable PEEP valve

Requirement		Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B	Flow approx L/min
40%	5	12 + 0	57
	7.5	12 + 0	52
	10	12 + 0	47
	12.5	12 + 0	42
50%	5	10 + 12	58
	7.5	10 + 12	54
	10	10 + 10	47
	12.5	10 + 10	42

Requirement		Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B	Flow approx L/min
60%	5	8 + 20	60
	7.5	8 + 20	55
	10	8 + 20	50
	12.5	9 + 16	42
100%	Feed with a minimum O ₂ flow of 42 L/min		

Performance table subject to changes. Always refer to product labelling.

Code	Description	Size	Box Qty.
CV200M	Ventumask CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Small / Medium	5
CV200L	Ventumask CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Large	5
CV200XL	Ventumask CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Extra Large	5
STS.15/30	Dual oxygen flow meter 15 + 30 L/min with O ₂ coupling		1



www.intersurgical.be/info/starmed

Video available at www.intersurgical.be

Make an enquiry

Dual flow meter manufactured by FLOW METER S.p.A.



Intersurgical Benelux B.V. Vluchttoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.be www.intersurgical.be



Manufacturer:
Intersurgical S.p.A.
Mirandola (MO) Italy

Please think before you print

Save energy and paper. If you must print this information sheet please print it double sided.

IS1.23_BENL • Issue 8 06.20

UK • Ireland • France • Germany • Spain • Portugal • Italy • Benelux • Sweden • Denmark • Lithuania • Russia • Czech Republic
Turkey • South Africa • China • Japan • Taiwan • Philippines • USA • Canada • Colombia • Australia