



## InterGuide tracheal tube introducers and InterForm intubating stylets

For use in the difficult airway




**Airway Management** ▪ Airway Accessories

## InterGuide tracheal tube introducer




InterGuide is a flexible tracheal tube introducer commonly known as a bougie. It allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations. InterGuide is single use, latex free and supplied sterile in easy-to-open individual packs.

### Atraumatic coudé tip

To reduce the potential for patient trauma



**Low-friction surface**  
For easy insertion and operation

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8070015	InterGuide tracheal tube introducer bougie	15FR	5.0mm	700mm	10 
8070010	InterGuide tracheal tube introducer bougie	10FR	3.3mm	700mm	10 
8070006	InterGuide tracheal tube introducer bougie	6FR	2.0mm	530mm	10 


Make an enquiry

## InterForm intubating stylet




InterForm is a malleable stylet, allowing the user to form the endotracheal tube into a suitable shape to ease insertion and give more control of the airway device. InterForm is single use, latex free and supplied sterile in easy-to-open individual packs.

### Aluminium core

Allows user to easily form suitable shapes



**Low-friction surface**  
For easy insertion and removal

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8080014	InterForm intubating stylet	14FR	4.7mm	340mm	10 
8080010	InterForm intubating stylet	10FR	3.3mm	340mm	10 
8080006	InterForm intubating stylet	6FR	2.0mm	275mm	10 

Make an enquiry

 Sterile



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland  
T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.

IS6.5 NL • Issue 4 01.18