## Intersurgical EcoLite<sup>™</sup> adult mask correct positioning



The elastic can be fitted either above or below the ears In the preferred under ear position it eliminates trauma to the ears



Code	Description	Box Qty.	
1136015	Intersurgical EcoLite, adult, medium concentration oxygen mask	55	
1135015	Intersurgical EcoLite, adult, medium concentration oxygen mask with tube, 2.1m	40	
1142015	Sentri Intersurgical EcoLite, adult, mask with CO2 monitoring line and tube, 2.1 m	30	
1143015	Sentri Intersurgical EcoLite, adult, mask with CO2 monitoring line, filter and tube, 2.1m	30	
1141015	Sentri Intersurgical EcoLite, adult, mask with tube, 2.1m	30	
1060085	Intersurgical EcoLite, adult, oxygen mask with 60% venturi valve, green	40	
1040085	Intersurgical EcoLite, adult, oxygen mask with 40% venturi valve, red	40	
1035085	Intersurgical EcoLite, adult, oxygen mask with 35% venturi valve, yellow	40	6
1031085	Intersurgical EcoLite, adult, oxygen mask with 31% venturi valve, orange	40	6
1028085	Intersurgical EcoLite, adult, oxygen mask with 28% venturi valve, white	40	6
1024085	Intersurgical EcoLite, adult, oxygen mask with 24% venturi valve, blue	40	
1107085	Intersurgical EcoLite, adult, venturi mask kit	25	
1181015	Intersurgical EcoLite, adult, high concentration oxygen mask with tube, 2.1m	24	
1188015	Intersurgical EcoLite, adult, aerosol mask	50	6
1453015	Cirrus2 nebuliser, adult, Intersurgical EcoLite mask kit with tube, 2.1m	30	
1463015	HOT Top2 nebuliser, adult, Intersurgical EcoLite mask kit with tube, 1.8m	30	6

Lower environmental impact product

ECOLITEPOS - Issue 1 10.18

INTERSURGICAL COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.co.uk www.intersurgical.co.uk