




## i-view™ video laryngoscope

Video laryngoscopy wherever and whenever you intubate



**Airway Management** ▪ Airway Accessories

## Introducing the i-view™ video laryngoscope

### One size ▪ Single use ▪ Fully disposable

i-view is the new, single use, fully disposable video laryngoscope from Intersurgical, providing the option of video laryngoscopy wherever you might need to intubate.

By incorporating a Macintosh blade, i-view can also be used for direct laryngoscopy and the technique for insertion is more familiar and instinctive than for devices with a hyper-angulated blade. Ready to use seconds after removing from the packaging, the ergonomic design ensures i-view is easy to use, and the integral LCD screen provides an optimal view in a variety of light conditions.



### i-view provides the opportunity of video laryngoscopy 'whenever and wherever you might need to intubate'

Where availability of a video laryngoscope may be limited due to the cost implications of purchasing reusable devices for multiple sites, i-view provides a cost effective solution, by combining all the advantages of a fully integrated video laryngoscope in a single use, disposable product. This makes i-view ideal for use in

- pre-hospital
- emergency medicine
- resuscitation
- anaesthesia
- difficult airway
- intensive care



Intensive care



Pre-hospital/resuscitation



Difficult airway



## Features & benefits

### Video laryngoscopy wherever and whenever you intubate

#### One size

no need to stock multiple sizes

#### Single use

reduces the risk of cross infection and avoids reprocessing costs

#### Fully disposable

no need to remove batteries prior to disposal

#### Integral LCD screen

provides an optimal view in a variety of light conditions

#### Ergonomic design

ensures i-view is comfortable and instinctive to use



#### On/Off switch

i-view is ready to use seconds after removing from the packaging

#### Camera and LED light

for optimal visualisation



#### Macintosh blade

i-view can also be used for direct as well as video laryngoscopy



Video available at [www.intersurgical.nl](http://www.intersurgical.nl)







i-view is supplied in boxes of 8 units. Each device is individually packed in a specially designed rigid shell sealed with a plastic film to ensure protection of the device, both in transit and in storage. The pack is easily opened by identifying the corner tab with the thumb symbol and pulling back the top plastic film until the product can be removed easily. The LOT number and expiry date of the device is located on the top edge above the LCD screen.



Code	Description	Box Qty.
8008000	i-view video laryngoscope	8

Make an enquiry



[www.intersurgical.nl/info/iview](http://www.intersurgical.nl/info/iview)



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland  
 T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.

IS6.15 NL ▪ Issue 1 03.18