



Silente™ adjustable venturi valve

Low operating noise



Oxygen and Aerosol Therapy ▪ Fixed Oxygen Concentration (High Flow)

Silente adjustable venturi valve

The Silente has been designed to operate at a very low noise level. The large diameter body provides space for clear, widely spaced calibration and incorporates a "pull-twist" action to prevent accidental adjustment.

The twin-jet design offers the user seven oxygen concentrations. The increased body diameter provides space for clear, widely spaced calibration and incorporates a 'pull-twist' action preventing accidental adjustment. The one-piece design eliminates the need to use more than one valve for all concentrations and there are no separate pieces to get lost or mislaid.

The concentration delivered is unaffected by the patient's breathing pattern provided that the optimum flow indicated on the valve is set on the oxygen flow meter.

Ease of use

The Intersurgical adjustable venturi valve is simple to use and has the following features and benefits:

- The appropriate flow range is selected to suit the desired concentration
- The outer ring with its pull-twist action selects the prescribed oxygen concentration
- A flow in L/min is printed above each concentration indicating the optimum flow meter setting.

Precision

The Silente delivers precise oxygen concentrations. This accuracy is independent of the oxygen flow meter setting and is achieved using just one valve.

Greener

The elimination of PVC from the masks has dramatically lowered their environmental impact

Soft face seals

Increased patient comfort

Elastic positioned under the ears

Eliminates trauma to top of ears

Pull-Twist controls

Prevent inadvertent adjustment

Make an enquiry

Incurved nose seal

Conforms to different nose shapes designed to prevent oxygen entering patient's eyes

Large, well spaced calibration settings

Easy to read, minimises errors





One piece design

No parts to change or lose

Extremely quiet

Much improved patient environment



Code	Description	Box Qty.
1137000	Silente, Eco, adult, adjustable venturi mask kit with tube, 1.8m	80 
1138015 NEW	Silente, Intersurgical EcoLite™, paediatric, adjustable venturi mask kit with tube, 1.8m	40 
1188015	Intersurgical EcoLite, adult, aerosol mask	50 
1198015	Intersurgical EcoLite, paediatric, aerosol mask	42 
1139000	Silente, adjustable venturi valve	50

 Lower environmental impact product



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

Please think before you print

Save energy and paper. If you must print this information sheet please print it double sided.

IS10.21 NL • Issue 4 04.18