




PEP Therapy Range Positive Expiratory Pressure



Oxygen and Aerosol Therapy • Respiratory Physiotherapy

The Intersurgical PEP therapy range

Intersurgical offer a range of products for the delivery of positive expiratory pressure therapy. The modular design of the product range enables clinicians to select from nine colour coded PEP resistors and two patient interfaces. The PEP T-piece is supplied with a mouthpiece that can easily be exchanged for a face mask if required. A pressure manometer can be connected to the system using the PEP manometer connector (1956020) designed to be fitted between the T-piece and PEP resistor.



Code	Description	Box Qty.
The PEP resistors are packaged in bags of 10 units		
8020060	PEP Resistor 6.0mm, purple, single patient use	100
8020050	PEP Resistor 5.0mm, brown, single patient use	100
8020045	PEP Resistor 4.5mm, red, single patient use	100
8020040	PEP Resistor 4.0mm, orange, single patient use	100
8020035	PEP Resistor 3.5mm, green, single patient use	100
8020030	PEP Resistor 3.0mm, blue, single patient use	100
8020025	PEP Resistor 2.5mm, yellow, single patient use	100
8020020	PEP Resistor 2.0mm, white, single patient use	100
8020015	PEP Resistor 1.5mm, black, single patient use	100
PEP T-piece, patient interfaces and connectors		
1956013	PEP T-piece, Inhalation port -22M/15F-22M with one-way valve and mouthpiece, 22F, single patient use	80
1956020	PEP manometer connector, 15M, 22M/15F, 6mm elbow with tube, 1.3m	40
7296001	ClearLite, anaesthetic face mask, size 6, extra large adult, red seal, no hook ring, 22F	20
7295001	ClearLite, anaesthetic face mask, size 5, large adult, orange seal, no hook ring, 22F	35
7294001	ClearLite, anaesthetic face mask, size 4, adult, green seal, no hook ring, 22F	35
7293001	ClearLite, anaesthetic face mask, size 3, small adult, yellow seal, no hook ring, 22F	35
7292001	ClearLite, anaesthetic face mask, size 2, paediatric, white seal, no hook ring, 22F	25

Make an enquiry

Watch the video

IS10.37 NL • Issue 3 08.21



Intersurgical Benelux B.V. Vluchttoord 28, 5406 XP, Uden, Nederland

T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016, ISO 14001:2015 and MDSAP

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.