



VT Sure Bag-Valve-Mask (BVM) resuscitators



Anaesthesia ▪ Resuscitation and Emergency Care




VT Sure Bag-Valve-Mask (BVM) resuscitators

In an emergency situation time is vital, there is little opportunity to check patient history and expensive re-usable equipment can often become contaminated. To avoid cross infection from blood, vomit and secretions we offer a complete range of single-use BVM resuscitators for use in resuscitation and emergency care.

The high quality, robust construction complies with international standard ISO 10651-4 and provides the reliability needed to withstand use in all clinical scenarios.

Supplied as a fully assembled single use product, the VT Sure BVM resuscitator is always ready for use in an emergency situation. Designed for ease of use, the VT Sure BVM resuscitator is available in a full range of sizes: Adult (1.5 L), small adult (1 L) and paediatric (550 ml).

The VT Sure range of BVMs are available with additional options such as PEEP valve, manometer or pressure relief valve to ensure safe and effective resuscitation.



BVM resuscitators are single-use and latex free.



1.5 L bag, size 5 mask

Code	Description	Box Qty.
7152111	VT Sure BVM resuscitator adult, 1.5 L bag with size 5 mask	8
7152121	VT Sure BVM resuscitator adult, 1.5 L bag with size 5 mask and manometer	8
7152131	VT Sure BVM resuscitator adult, 1.5 L bag with size 5 mask, manometer and adjustable PEEP valve	8

1 L bag, size 4 mask

Code	Description	Box Qty.
7153111	VT Sure BVM resuscitator small adult/paediatric, 1 L bag with pressure relief valve 40 cm/H ₂ O and size 4 mask	8
7153121	VT Sure BVM resuscitator small adult/paediatric, 1 L bag with pressure relief valve 40 cm/H ₂ O, size 4 mask and manometer	8
7153131 	VT Sure BVM resuscitator small adult/paediatric, 1 L bag with pressure relief valve 40 cm/H ₂ O, size 4 mask, manometer and adjustable PEEP valve	8

550 ml bag, size 3 mask

Code	Description	Box Qty.
7151111	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O and size 3 mask	8
7151121	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O size 3 mask and manometer	8
7151131	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O size 3 mask, manometer and adjustable PEEP valve	8

550 ml bag, size 1 mask

Code	Description	Box Qty.
7154111	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O and size 1 mask	8
7154121	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O, size 1 mask and manometer	8
7154131	VT Sure BVM resuscitator paediatric, 550 ml bag with pressure relief valve 40 cm/H ₂ O, size 1 mask, manometer and adjustable PEEP valve	8

 Illustrated

VT Sure Bag-Valve-Mask is manufactured by Pulmodyne Inc and distributed by Intersurgical Ltd.



Accessories



7162111



7162112



7163000



1644000

Code	Description	Box Qty.
7162111	Dial pressure manometer	10
7162112	Angled pressure manometer	10
7163000	BVM mask strap	10
1644000	Clear-Guard™ Midi low volume breathing filter with luer port	100