

# STARMed CaStar

Hood for CPAP therapy



# STARMed CaStar – CPAP Hood

## StarMed CaStar

The CaStar is a comfortable, versatile and lightweight CPAP patient interface designed to offer an alternative to traditional methods of CPAP delivery. The CaStar is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in both semi-recumbent and supine positions providing an ideal option for CPAP therapy.

The CaStar may be used for the treatment of a number of respiratory conditions including:

Pre and post extubation respiratory support, COPD exacerbation, thoracic trauma, asthma, post operative hypoxaemia, cardiogenic and non cardiogenic pulmonary oedema.

## Technical Features

- 22M/22F connectors compliant to ISO standard
- Fitted with integrated pressure gauge to check the presence of pressure within the hood
- Internal volume 16 litres – weight 340g
- Comfortable and atraumatic underarm fastening straps
- Other fastening accessories are available if required
- Two sealed access ports for probes or catheters of 3.5-7.0mm Ø
- Patient access port with screw cap. Fitted with bidirectional anti-suffocation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated



## Features and benefits



CE certified product

### Transparent

For ease of patient monitoring

### Integrated pressure manometer

### Patient access port with anti-suffocation valve

### Comfortable underarm fastening

CaStar CPAP hood with anti-suffocation valve, patient access port, ear plugs and measuring tape

Code	Size	Neck Circumference	Box qty.
<b>CP210XXL/2W</b>	Extra Extra Large	≥50cm	5
<b>CP210XL/2W</b>	Extra Large	45-52cm	5
<b>CP210L/2W</b>	Large	40-47cm	5
<b>CP210M/2W</b>	Medium	34-41cm	5
<b>CP210S/2W</b>	Small	27-34cm	5
<b>CP210XS/2W</b>	Extra Small	17-27cm	5

Single packed. Minimum purchasing package: 5 units. Custom-packs are supplied on specific request of the customer.

IS1.17\_UK • Issue 3\_05.13



Crane House, Molly Millars Lane,  
Wokingham, Berkshire, RG41 2RZ  
T. +44 (0)118 9656 300  
F. +44 (0)118 9656 356

UK  
info@intersurgical.co.uk  
Italy  
info@intersurgical.it

España  
info@intersurgical-es.com  
Portugal  
info@intersurgical.pt

Lietuva  
info@intersurgical.lt  
Россия  
info@intersurgical.ru

France  
info@intersurgical.fr  
South Africa  
info@intersurgical.co.za

Japan  
info@intersurgical.co.jp  
Nederland  
info@intersurgical.nl



Manufactured by: Intersurgical S.p.A.