



Fabian HFO and fabian evolution neonatal ventilator range


Breathing systems suitable for use with Acutronic fabian ventilators



Critical Care ▪ Breathing Systems

Introducing a comprehensive range of breathing systems suitable for use with the Acutronic fabian HFO and the fabian evolution ventilators, offering conventional and high frequency ventilation choices. Using both Smoothbore and Flextube™ tubing, the range features dual and single heated wire systems, infant nasal CPAP system as well as a transport option.

We have included an nFlow heated wire breathing system and an nFlow adaptor kit to convert any of our featured breathing systems to fit the fabian HFO and fabian evolution ventilators. With a new detachable monitoring line allowing for out-of-incubator care we now offer a complete care solution for infant nasal CPAP.



	Code	Description	Box Qty.
👁️	4651000	HFO adaptor kit to modify fabian breathing systems to HFO ventilation	20
	4650000	10mm Flextube transport breathing system with swivel Y-piece and monitor line for fabian ventilators, 1.2m	15
	4652000	10mm Flextube transport breathing system with swivel Y-piece and monitor line for fabian ventilators, 1.6m	10
	4760000	nFlow infant nasal CPAP adapter kit with detachable monitoring line for fabian ventilators, 1.6m	15
	4761000	nFlow infant nasal CPAP Flextube breathing system with heated wire for fabian ventilators, 1.6m	10
	4071000	10mm Flextube single heated wire breathing system with swivel Y-piece, water trap and monitoring line for fabian ventilators, 1.6m	12
👁️	4072000	10mm Flextube dual heated wire breathing system with swivel Y-piece and monitoring line for fabian ventilators, 1.6m	12
	5071000	10mm Smoothbore single heated wire breathing system with swivel Y-piece, water trap and monitoring line for fabian ventilators, 1.6m	20
	5072000	10mm Smoothbore dual heated wire breathing system with swivel Y-piece and monitoring line for fabian ventilators, 1.6m	20

IS1.38 NL • Issue 1 01.14