

7000 Series Respiratory Humidifier



Introducing the new 7000 Series Respiratory Humidifier from Intersurgical

Intersurgical is pleased to introduce the new 7000 Series Respiratory Humidifier to the range. This will provide a further customer option for active humidification, allowing us to offer a complete solution in this clinical area, alongside our extensive range of breathing systems and accessories for adults, paediatrics and neonates.

The 7000 Series Respiratory Humidifier offers a number of key benefits:

- Simple intuitive device whilst providing advanced options.
- Includes a Manual mode option, allowing users the ability to pre-set airway and chamber temperature.
- The user has the choice of which temperature is displayed as default option.
- An Alarm indicator panel with icon display of alarm location and priority.

The 7000 Series Respiratory Humidifier is compatible with Intersurgical humidification chambers and heated wire breathing systems offering a complete solution

Features and Benefits

Locking Button

One touch mounting allows for easy insertion and removal of the humidification chamber

Configurable Display

Choice of default either chamber or patient end airway temperature

Mode Selection

IPPV, NIV and a Manual mode which allows the user to set the desired temperatures

Alarm Indicator Panel

Highlights alarm location and priority





Ordering Information




7000 Series Respiratory Humidifier

		
7000000	Respiratory Humidifier (UK)	Box qty. 1
7000001	Respiratory Humidifier (JPN)	Box qty. 1
7000002	Respiratory Humidifier (EU1)	Box qty. 1

Humidification Chambers

<div>  <div>2310000</div> </div> <div>  <div>2330000</div> </div>		
2310000	Auto-fill humidification chamber	Box qty. 30
2330000	Low volume manual fill humidification chamber	Box qty. 30

Accessories

		
5601000	Electrical adaptor lead (Dual heated wire system)	Box qty. 1
		
5600000	Electrical adaptor lead (Single heated wire system)	Box qty. 1
		
7700000	Dual temperature probe assembly	Box qty. 1

Breathing systems for Active humidification with heated wire

system illustrated

2026000	22mm Flextube, Dual Heated wire breathing system, 0.5 limb, 1.6m	Box qty. 10
2026310	22mm Flextube, Dual Heated wire breathing system, with autofill chamber, 0.5m limb, 1.6m	Box qty. 7

2025310	22mm Flextube, Single Heated wire breathing system, with autofill chamber, water trap and 0.5m limb, 1.6m	Box qty. 7
----------------	-----------------------------------------------------------------------------------------------------------	------------

We have a wide range of breathing systems for active humidification for adults, paediatrics and neonates, please refer to:
www.intersurgical.com/catalogue - for full range Product Catalogue (Edition 12)
www.intersurgical.com/neonatal-catalogue - for Neonatal Catalogue (Edition 7)

IS1.29 Issue 3 04.13