





InterGuide tracheal tube introducers and InterForm intubating stylets

For use in the difficult airway




Airway Management ▪ Airway Accessories

InterGuide tracheal tube introducer

InterGuide is a flexible tracheal tube introducer commonly known as a bougie. It allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations. InterGuide is single use, latex free and supplied sterile in easy-to-open individual packs.

Atraumatic coude tip

To reduce the potential for patient trauma



Low-friction surface
For easy insertion and operation

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8070015	InterGuide tracheal tube introducer bougie	15FR	5.0mm	700mm	10
8070010	InterGuide tracheal tube introducer bougie	10FR	3.3mm	700mm	10
8070006	InterGuide tracheal tube introducer bougie	6FR	2.0mm	530mm	10

InterForm intubation stylet

InterForm is a malleable stylet, allowing the user to form the endotracheal tube into a suitable shape to ease insertion and give more control of the airway device. InterForm is single use, latex free and supplied sterile in easy-to-open individual packs.


Aluminium core

Allows user to easily form suitable shapes



Low-friction surface
For easy insertion and removal

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8080014	InterForm intubation stylet	14FR	4.7mm	340mm	10
8080010	InterForm intubation stylet	10FR	3.3mm	340mm	10
8080006	InterForm intubation stylet	6FR	2.0mm	275mm	10

 All products are supplied sterile

IS6.5 • Issue 2 03.14



Intersurgical Nederland B.V, Vluchtoord 28, 5406 XP, Uden, Nederland
T: 0413 243 860 F: 0413 243 869 info@intersurgical.nl www.intersurgical.nl

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy
Netherlands ▪ Sweden ▪ Lithuania ▪ Russia ▪ Czech Republic
South Africa ▪ China ▪ Japan ▪ Taiwan ▪ Philippines
USA ▪ Colombia ▪ Australia



Please think before you print
Save energy and paper.
If you must print this e-mail
please print it double sided

Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485