



One-piece Guedel airways

Safer by design



Airway Management ▪ Airway Devices

One-piece Guedel airways

Safer by design

By using the latest manufacturing technology Intersurgical has produced the first major breakthrough in Guedel airway design for thirty years.


Combining two materials with very different physical characteristics into one moulding, our One-piece design eliminates the separate bite block found in all conventional Guedel airways.

Loose or detached bite blocks constitute a real threat to patient safety and have been the subject of several hazard notices in the past. Our One-piece design overcomes all of the associated risks.

As well as improved safety, the One-piece also incorporates a soft tip reducing the risk of trauma to the patient. We also offer a sterile option, providing a further choice.



▶ Video available at www.intersurgical.nl



Rigid body

To maintain a patent airway

Soft tip

Reduces risk of tissue damage

One-piece design
With no separate inner bite block – safer by design

Nine colour-coded sizes

To meet all patient requirements



Softer second material
Reduces the potential for dental damage

Clear colour-coding
For easy identification

Code	Description	Size	ISO	Colour	Box Qty.
1115120 (S*)	One-piece Guedel airway	5	12.0	Purple	70
1114100 (S*)	One-piece Guedel airway	4	10.0	Red	80
1113090 (S*)	One-piece Guedel airway	3	9.0	Orange	80
1112080 (S*)	One-piece Guedel airway	2	8.0	Green	100
1111570 (S*)	One-piece Guedel airway	1.5	7.0	Yellow	100
1111065 (S*)	One-piece Guedel airway	1	6.5	White	20
1110055 (S*)	One-piece Guedel airway	0	5.5	Grey	25
1100050 (S*)	One-piece Guedel airway	00	5.0	Blue	25
1000035 (S*)	One-piece Guedel airway	000	3.5	Pink	25

Make an enquiry

Lower environmental impact product

Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland
T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.