



nFlow™

Infant nasal CPAP system



Critical Care ▪ Breathing Systems



nFlow™

Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

A choice for infant CPAP

The Intersurgical nFlow system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.

nFlow is compatible with flow drivers designed for use with gas flip technology CPAP, providing slightly higher system pressure associated with this treatment.

The nFlow system

The nFlow system is supplied assembled with generator, heated wire limb and 10mm adjustable Superset style exhaust tube for efficient gas removal. The complete nFlow breathing system is phthalate free including the small bore tube used for monitoring and gas delivery.

A choice of patient interfaces and bonnets



Patient interfaces and bonnets are supplied separately to avoid waste. Use the size guide provided to choose the correct option for your patient.



Make an enquiry

Choice made easy

The nFlow system includes two sizing guides making it really simple to choose the correct fitting interface and bonnet for the comfort of your patient.



12 bonnets
Available separately in twelve colour-coded sizes. Each bonnet includes three Velcro® fasteners and secure attachment points.

Patient interface
Masks and prongs are available separately in four size options. Made from soft, comfortable silicone to ensure optimal patient fit with minimal wastage.





We've made ordering easy. To make a complete system you can choose one product from group A and B, and choose a patient interface from C and/or D. At any time you can order separate items according to your needs.



| Code | Description | Box Qty. |
|------|-------------|----------|
|------|-------------|----------|

A nFlow breathing systems

| | | | |
|--|---------|--|----|
| | 4700000 | nFlow infant CPAP heated wire breathing system, 1.6m | 12 |
| | 4700001 | nFlow infant CPAP heated wire breathing system, 1.8m | 10 |
| | 4700016 | nFlow infant CPAP heated wire breathing system with detachable monitoring line, 1.6m | 10 |
| | 4708000 | nFlow silencer kit | 50 |
| | 4700011 | nFlow replacement bonnet ties | 60 |

B nFlow bonnets Size

| | | | | |
|--|---------|------------------------------------|-----------------------|----|
| | 4707009 | nFlow, infant bonnet, navy | size 9 (40 - 42 cm) | 10 |
| | 4707008 | nFlow, infant bonnet, turquoise | size 8 (38 - 40 cm) | 10 |
| | 4707007 | nFlow, infant bonnet, orange | size 7 (36 - 38 cm) | 10 |
| | 4707006 | nFlow, infant bonnet, red | size 6 (34 - 36 cm) | 10 |
| | 4707005 | nFlow, infant bonnet, green | size 5 (32 - 34 cm) | 20 |
| | 4707004 | nFlow, infant bonnet, light orange | size 4 (30 - 32 cm) | 20 |
| | 4707003 | nFlow, infant bonnet, blue | size 3 (28 - 30 cm) | 20 |
| | 4707002 | nFlow, infant bonnet, yellow | size 2 (26 - 28 cm) | 20 |
| | 4707001 | nFlow, infant bonnet, dark grey | size 1 (24 - 26 cm) | 20 |
| | 4707110 | nFlow, infant bonnet, pink | size 0 (22 - 24 cm) | 20 |
| | 4707100 | nFlow, infant bonnet, grey | size 00 (20 - 22 cm) | 20 |
| | 4707000 | nFlow, infant bonnet, white | size 000 (18 - 20 cm) | 20 |

C nFlow nasal masks Size

| | | | | |
|--|---------|------------------|-------------|----|
| | 4716000 | nFlow nasal mask | extra large | 30 |
| | 4706000 | nFlow nasal mask | large | 30 |
| | 4705000 | nFlow nasal mask | medium | 30 |
| | 4704000 | nFlow nasal mask | small | 30 |

D nFlow nasal prongs Size

| | | | | |
|--|---------|-------------------|-------------|----|
| | 4703000 | nFlow nasal prong | large | 30 |
| | 4702000 | nFlow nasal prong | medium | 30 |
| | 4701000 | nFlow nasal prong | small | 30 |
| | 4711000 | nFlow nasal prong | extra small | 30 |



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland
 T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



The manufacturer Intersurgical Ltd is certified to ISO 9001:2008, ISO 13485:2003 and ISO 14001:2015

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.

IS6.16 NL • Issue 7 01.19