



STARMed CaStar Up

Openable hood for CPAP therapy



Critical Care ▪ Patient Interfaces ▪ **STARMed**

STARMed CaStar Up openable CPAP hood

StarMed CaStar Up

Openable CaStar Up is a comfortable, versatile and lightweight CPAP patient interface designed to offer an alternative to traditional methods of CPAP delivery. The CaStar Up is well tolerated for extended periods and can reduce the risks associated with both endotracheal intubation and face mask use. It can be used in semi-recumbent, supine and prone positions providing an ideal option for CPAP therapy.

Advantages

- Allows rapid access to the patient during nursing procedures or in case of emergency
- One operator can easily remove the upper part of the hood and have full access to the patient
- Probes and catheters can remain in place, on the patient, connected to the lower part of the hood

Technical Features

- Manufactured from phthalate-free materials
- 22M/22F connectors compliant to ISO standard
- Fitted with integrated pressure gauge to check the presence of pressure within the hood
- Built in inflatable neck cushion for added patient comfort
- Comfortable and atraumatic underarm fastening straps
- Other fastening accessories are available if required
- Two versatile sealed access ports suitable for single or multiple lumen probes or catheters
- New design patient access port. Fitted with bi-directional anti-asphyxiation valve which opens automatically in the event of pressure failure
- Airtight sealing of the system is achieved by an elastic membrane making contact with the upper part of the torso, once the hood is inflated

Read more evidence on www.intersurgical.nl

CaStar Up openable CPAP hood with inflatable neck cushion, bi-directional anti-asphyxiation valve and patient access port*

Code	Size	Neck circumference	Box Qty.
CA012XXL/2	Extra Extra Large	≥50cm	5
CA012XL/2	Extra Large	45-52cm	5
CA012L/2	Large	40-47cm	5
CA012M/2	Medium	34-41cm	5
CA012S/2	Small	27-34cm	5
CA012XS/2	Extra Small	17-27cm	5



www.intersurgical.nl/info/starmed

Video available at www.intersurgical.nl

Make an enquiry

*PEEP valve and tubing not included and have to be ordered separately. Please refer to the website for the range of accessories.



Intersurgical Benelux B.V. Vluchtoord 28, 5406 XP, Uden, Nederland
 T: +31 (0)413 243 860 F: +31 (0)413 243 869 info@intersurgical.nl www.intersurgical.nl



Manufacturer:
Intersurgical S.p.A.
Mirandola (MO) Italy

Please think before you print
Save energy and paper. If you must print this information sheet please print it double sided.

IS1.31_NL • Issue 7 06.20